CLEAR ALIGNER

22

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All cases were carried out by Dr. Pablo Echarri. All the laboratory work was carried out by Ladent laboratory, and the Scheu Dental Technology materials were used.

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© 2015 Centro de Ortodoncia y ATM, Ladent, SL
All rights reserved

This book or any part thereof may not be reproduced, stored in retrieval system or transmitted in any form or by any means electronic, mechanical, photocopying or otherwise, without prior written permission of the publishers.

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CA Clear Aligner anchorage and retention

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Figure 1. In CA Clear Aligner technique, the most frequently treated malpositions are the incisors and canines malpositions, and this is why we distinguish the treatment zone (A), incisors & canines zone, and the retention and anchorage zone (B), bicuspid & molars zone. The anchorage capacity is essential for any orthodontic appliance, because it defines the treatment possibilities.

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
Figure 2. Distal zone of the last erupted molar is very important for the CA Clear Aligner retention. Because of this, it is very important that this zone is clearly defined in impressions and in models sent to the laboratory. When the third molars erupt, they reduce the retention of aligners.

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
Figure 3. Incomplete eruption of the last molar in the arch, and especially of its distal surface, reduces the retention capacity of aligners.

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Figure 4. The anatomical shape of molar crowns is very important. Cone-shaped crowns offer much less retention than more retentive-shaped crowns.

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Figure 5. Molars with alterations in shape due to fractures or caries, or with furca lesions offer less retention.

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Figure 6. Molars with badly adjusted prosthetic crowns offer less retention than when the prosthetic crowns perfectly fit.

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Figure 7. Molars with crowns of reduced height offer much less retention than molars with high crowns.

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Figure 8. The patients with absent molars offer less retention and anchorage than the patients with all molars.

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Figure 9. If there are spaces between molars, the retention capacity of aligners is reduced.

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Figure 10. The patients with molars with low cusps, or cusps presenting abrasion offer less retention of aligners due to the reduction of molar height. Furthermore, the lack of intercuspsation reduces the anchorage offered by these molars.

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Figure 11. There are two ways to increase the retention of aligners. The first is to create an extra retention in the crown, and for this purpose CA Power Grips nº 4 FALA (Force Application Labial) or nº 5 FALI (Force Application Lingual) are used. The one which fits better to the tooth will be used.

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Figure 12. Shape of a composite retention carried out with FALA. For clinical procedure of FALA application, please see CA Newsletter nº8.

Figure 13. Shape of a composite retention carried out with FALI. For clinical procedure of FALI application, please see CA Newsletter nº8.

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Figure 16. Another way of increasing the CA Clear Aligner retention is to use CA Plier Tip 1, as shown in the drawing.

Figure 17. CA Plier Tip 1, used to increase the "cold" retention (it is not necessary to heat up the plier).

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Figure 18. CA Plier Tip 1, active part of the plier.

Figure 19. How to use the CA Plier Tip 1 to increase the retention of CA Clear Aligner.

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
Figure 20. Additional retention in a CA Clear Aligner carried out with the plier between the contact point and interdental papilla.

Figure 21. In sagittal or oblique-sagittal expansion treatments the aligner should apply the reciprocal force towards distal over the posterior teeth.

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
Figure 22. In space closure treatments with retroclination, the aligner will tend to carry out the reciprocal force towards mesial over the posterior teeth.

Figure 23. In anterior teeth intrusion treatments there are no retention problems in the anterior part of the aligner. In anterior deep bite treatments, it is important to increase the retention in the posterior part of the aligner.

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
Figure 24. In anterior open bite treatments in which the anterior teeth extrusion is indicated, the retention of anterior part of the aligner is not sufficient, and it is necessary to indicate the use of CA Power Grip nº1 in anterior teeth.

Figure 25. Ca Clear Aligner for extrusion of anterior teeth with CA Power Grips in incisors.

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SITUATION	SOLUTION	
	CA Power Grip 4 or 5	CA Plier Tip 1
Eruption of the third molars	X	X
Eruption of the distal part of 2nd molars	X	X
Cone-shaped molar crowns	X	X
Caries or furcae	Treatment of caries or furcae	
Badly adjusted prosthetic crowns	Replacement of prosthetic crowns or CA Power Grips 4 or 5	
Short crowns of molars	X	X
Absent molars	X	
Spacing between molars	X	
Low cusps, or cusps with abrasion	X	X