CLEAR ALIGNER

ACTUALIZANDO CA[®] CLEAR ALIGNER
CON LOS TORNILLOS CA[®] VECTOR[®]

UPDATING CA[®] CLEAR ALIGNER
WITH CA[®] VECTOR[®] SCREWS

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All cases were carried out by Dr. Pablo Echarri. All the laboratory work was carried out by Ladent laboratory, and the Scheu Dental Technology materials were used.

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© 2017 Centro de Ortodoncia y ATM, Ladent, SL All rights reserved

This book or any part thereof may not be reproduced, stored in retrieval system or transmitted in any form or by any means electronic, mechanical, photocopying or otherwise, without prior written permission of the publishers.

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VECTOR® EXPANSION SCREWS FOR CA® CLEAR ALIGNER

TRANSVERSAL EXPANSION	ANTERIOR	UPPER	40 - 45
		LOWER	40
	POSTERIOR	UPPER	50 - 140 - 160
		LOWER	
LATERAL OBLIQUE	UPPER	(x2) 30 + CA Arch	
ONE-TOOTH EXPANSION	UPPER		PISTON
	LOWER		
DISTALIZATION	UPPER		350 DISTALIZER KIT
	LOWER		
CLOSING-SPACES	UPPER		TRACTIVE 350
	LOWER		

Fig. 1. Scheme of VECTOR® screws used in CA® CLEAR ALIGNER technique.

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ACTIVACIÓN

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INDICATION

Upper arch and lower arch transversal expansion.

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Upper arch transversal expansion.

POSITION

At the level of bicuspids in the middle line of the palate.

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Fig 8. Scheme of VECTOR® 30 screw with CA® Arch


Fig 9. Tornillo VECTOR® 30.

Fig 9. VECTOR® 30 screw.

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ACTIVACIÓN

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INDICATION

Lateral-oblique expansion in the upper arch

POSITION

The expansion Kit includes two one-way screws adapted in the stainless-steel arch.

OBJETIVE

Greater expansion in anterior zone.

ACTIVATION

1/4 turn of the screw each 2 to 3 days.

TORNILLO PISTON


Fig 10. Esquema del tornillo PISTON

Fig 11. Tornillo Piston

INDICACIÓN

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POSICIÓN

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OBJETIVO

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ACTIVACIÓN

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PISTON SCREW


Fig 10. Scheme of Piston Screw


Fig 11. Piston screw.

INDICATION

One-tooth expansion.

POSITION

In lingual zone of the tooth.

OBJETIVE

One-tooth expansion.

ACTIVATION

1/4 turn of the screw each 2 to 3 days.

CLINICAL PAPER

CA[®] DISTALIZER KIT
+
VECTOR[®] 350


Fig 12. Esquema del tornillo CA[®] Distalizer Kit + VECTOR[®] 350

CLEAR ALIGNER

CA DISTALIZER KIT
+
VECTOR[®] 350


Fig 12. Scheme of CA[®] Distalizer Kit + VECTOR[®] 350

Fig 13. CA[®] Distalizer Kit + VECTOR[®] 350

Fig 13. CA[®] Distalizer Kit + VECTOR[®] 350

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ACTIVACIÓN

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INDICATION

Upper arch and lower arch molar distalization

POSITION

One-way screw in lingual zone and tube, sectional arch and open coil-spring in buccal zone.

OBJECTIVE

Upper arch and lower arch molar distalization

ACTIVATION

1/4 turn of the screw each 2 to 3 days.

**TORNILLO TRACTIVE VECTOR®
350**


Fig 14. Esquema del tornillo Tractive VECTOR® 350

**TRACTIVE VECTOR® 350
SCREW**


Fig 14. Scheme of Tractive VECTOR® 350 screw

Fig 15. Tornillo Tractive VECTOR® 350.

Fig 15. Tractive VECTOR® 350 screw.

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ACTIVACIÓN

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INDICATION

Mesialization of upper and lower molars.

POSITION

One-way screw in lingual zone and tube, sectional arch and open coil-spring in buccal zone.

OBJETIVE

Mesialization of upper and lower molars.

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USAGE PROTOCOL OF VECTOR® SCREW


Fig. 16. Squeme of the Usage Protocol of VECTOR® screws