CLEAR ALIGNER

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USE PROTOCOL FOR
CA VECTOR SCREW WITH
CA® CLEAR ALIGNER DIGITAL

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All cases were carried out by Dr. Pablo Echarri. All the laboratory work was carried out by Ladent laboratory, and the Scheu Dental Technology materials were used.

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© 2018 Centro de Ortodoncia y ATM, Ladent, SL All rights reserved

This book or any part thereof may not be reproduced, stored in retrieval system or transmitted in any form or by any means electronic, mechanical, photocopying or otherwise, without prior written permission of the publishers.

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VECTOR® EXPANSION SCREWS FOR CA® CLEAR ALIGNER

TRANSVERSAL EXPANSION	ANTERIOR	UPPER	40 - 45
		LOWER	40
	POSTERIOR	UPPER	50 - 140 - 160
		LOWER	
	LATERAL OBLIQUE	UPPER	(x2) 30 + CA Arch

ONE-TOOTH EXPANSION	UPPER		PISTON
	LOWER		

DISTALIZATION	UPPER		350 DISTALIZER KIT
	LOWER		

CLOSING-SPACES	UPPER		TRACTIVE 350
	LOWER		

Fig. 3. Scheme of VECTOR® screws used in CA® CLEAR ALIGNER technique.

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ACTIVACIÓN

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INDICATION

Upper arch and lower arch transversal expansion.

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Upper arch transversal expansion.

POSITION

At the level of bicuspid in the middle line of the palate.

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+

VECTOR® 350


Fig 8. Scheme of CA® Distalizer Kit + VECTOR® 350

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Fig 9. CA® Distalizer Kit + VECTOR® 350

First distalize one site and then distalize the other side depending on the anchorage

Maximum apex movement is 3 mm (MD direction)

Fig 9. CA® Distalizer Kit + VECTOR® 350

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ACTIVACIÓN

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INDICATION

Upper arch and lower arch molar distalization

POSITION

One-way screw in lingual zone and tube, sectional arch and open coil-spring in buccal zone.

OBJECTIVE

Upper arch and lower arch molar distalization

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VECTOR® 30 SCREW WITH CA® ARCH


Fig 10. Scheme of VECTOR® 30 screw with CA® Arch



Fig 11. Tornillo VECTOR® 30.

Fig 11. VECTOR® 30 screw.

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ACTIVACIÓN

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INDICATION

Lateral-oblique expansion in the upper arch

POSITION

The expansion Kit includes two one-way screws adapted in the stainless-steel arch.

OBJECTIVE

Greater expansion in anterior zone.

ACTIVATION

1/4 turn of the screw each 2 to 3 days.

TORNILLO TRACTIVE VECTOR® 350


Fig 12. Esquema del tornillo Tractive VECTOR® 350

TRACTIVE VECTOR® 350 SCREW


Fig 12. Scheme of Tractive VECTOR® 350 screw

Fig 13. Tornillo Tractive VECTOR® 350.

Fig 13. Tractive VECTOR® 350 screw.

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ACTIVACIÓN

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INDICATION

Mesialization of upper and lower molars.

POSITION

One-way screw in lingual zone and tube, sectional arch and open coil-spring in buccal zone.

OBJECTIVE

Mesialization of upper and lower molars.

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USAGE PROTOCOL OF VECTOR® SCREW


Fig. 18. Squeme of the Usage Protocol of VECTOR® screws

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USE OF CA VECTOR SCREW WITH CA® CLEAR ALIGNER DIGITAL


Fig. 19. Registro inicial.

Fig. 19. Initial register.

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Fig. 20. Register after the use of CA Vector Screw for lateral-oblique expansion.

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If any CA Vector Screw is required, firstly is indicated the use of CA Vector Screw and, after the activation, it's necessary to send to the laboratory new registers for realization of final set-up.

Thus the technician is able to plan the treatment and the number of phases more accurate due to the use of real spaces.