CLEAR ALIGNER

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All cases were carried out by Dr. Pablo Echarri. All the laboratory work was carried out by Ladent laboratory, and the Scheu Dental Technology materials were used.

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©2016 Centro de Ortodoncia y ATM, Ladent, SL
All rights reserved

This book or any part thereof may not be reproduced, stored in retrieval system or transmitted in any form or by any means electronic, mechanical, photocopying or otherwise, without prior written permission of the publishers.

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The inter-proximal contact point is very important to keep the alignment of dental arches. If there is some crowding or rotations, the inter-proximal relationship is affected and the most likely is that the crowding and/or the rotations are going to increase.

Also the inter-proximal contact points facilitate the transmission of occlusal forces to the dental arch. Each tooth supports the neighboring teeth and when this balance is broken, modifications take place in the periodontal tissues.

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Fig. 1. Interdental fibres in the base of the papilla. These contact points keep the health of the dento-dental and supra-septal fibres.

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Immediatly inside the fibres which contribute to the approachment of the teeth, positioned in the base of interdental papillae, a peridontal insertion, very important for teeth preservation, can be found.

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For this reason it is very important that the orthodontists know the inter-proximal contact points from the anatomical point of view.

The contact point is usually in the occlusal third of the mesial or distal surfaces of all the teeth, but it also depends on the dental shape.

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Fig. 2. Shapes of the teeth: A. Rectangular, B. Triangular, C. "Barrel"-shaped

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In the rectangular teeth, the contact point is bigger and more gingival; in the triangular ones, it is smaller and in a more incisal position; and in the "barrel-shaped" teeth, the contact point is in the middle third in the incisal-gingival sense.

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Fig. 3. Position of the inter-proximal contact point in different teeth.

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In the anterior teeth, incisors and cuspids, the contact point is in the middle third in labio-lingual sense but in the posterior teeth, normally it is in the labial third.

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Fig. 4. Upper and lower arches and inter-proximal contact point.

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Fig. 5. The contact point can increase its surface and transform into an ovoid inter-proximal contact surface in all cases, but more vertical in anterior teeth and more horizontal in posterior teeth.

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According to Begg, the spontaneous mesial movement of the teeth is to compensate the mesial and distal surfaces abrasion during the mastication of hard food.

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These spaces are very important during the mastication because the movement of the food is sliding over these surfaces.

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Fig. 7. Movement of the food during the mastication: A - Mesial view: the food slides over the labial and lingual surfaces of the teeth (white arrows) and over the papilla (black arrows). B - Occlusal view: the food slides over the occlusal sulcus to the occlusal surfaces of the teeth and the labial and lingual sulcus. C - Labial view: circulation of the food. (Taken and modified from Aprile and Figún).

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The functions of interproximal contact point are to:

- Maintain the position of the teeth in the dental arch, avoiding crowding and rotations.
- Protect the dental papilla and the dento-dental and supra-septal fibres.
- Avoid the food retention.
- Aesthetics.

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Fig. 8. The length of the papilla is normally 4.5 mm. (Tarnow).

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Fig. 9. Relationship between the distance of the contact point and the bone crest with the presence and absence of the papilla.

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Therefore, if inter-proximal contact point is established 5 mm from the bone crest, the papilla is present in a 100% of cases. If the contact point is 6 mm from the bone crest, the papilla is present in only 56% of cases and if the same distance is increased up to 7 mm or more, the papilla is present in only 27% of cases (Tarnow).

Stripping and approximation allows the correction of gingival "black triangles" because they fix the contact point to 5 mm from the bone crest.

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Fig. 10. Correction of the gingival "black triangles" depending on the dento-alveolar discrepancy and the Bolton Index.

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Stripping and Bolton Index

The Bolton Index is the relation of the size in upper and lower teeth. The Bolton Index "12" takes into account the measures of the teeth from first molar to first molar and the Bolton Index "6" takes into account the measures of the teeth from canine to canine measurements.


Fig. 11. Índice de Bolton "12".

Fig. 11. Bolton Index "12".


Fig. 12. Índice de Bolton "6".

Fig. 12. Bolton Index "6".

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The following situations can be found:

If Bolton Index "12" is normal, it means that the case after the orthodontic treatment can be finished in molar Class I, and if Bolton Index "6" is normal the case can be finished in canine Class I. This would be the ideal situation and in another words means that there is no size discrepancy between the teeth of the upper arch and the teeth of the lower arch considering all the teeth from the upper first molar to the upper first molar of the other side for the Bolton Index "12" or having into account only the teeth from canine to canine for the Bolton Index "6".

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Fig. 13. Case without Bolton discrepancy.

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If there is Bolton index "12" excess of, for example, 4 mm, and there is also Bolton Index "6" excess of 4 mm too, i.e. both Bolton Indexes have equal excess, a canine-to-canine stripping in maxilla should be carried out, because the Bolton Index "6" is included into the Bolton Index "12", which means that the excess is only in the anterior sector.

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Fig. 14. Case with the maxillary excess of Bolton Index "12" and "6".

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If there is Bolton Index "12" excess of, for example, 4 mm, and Bolton Index "6" is normal, in this situation the upper molars and bicuspid should be stripped, because the excess is located in the posterior sector.

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Fig. 15. Case with excess of maxilla Bolton Index "12" and without Bolton Index "6" discrepancy.

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If there is Bolton Index "12" excess of, for example, 4 mm, and there is also Bolton Index "6" excess of, for example, 2 mm, i.e. Bolton Index "12" excess is higher than Bolton Index "6" excess, total stripping should be carried out. In this case there is an excess in both anterior and posterior sectors.

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Fig. 16. Case with maxilla Bolton Index "12" excess higher than Bolton Index "6" excess.

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The same reasoning is applied to Bolton Index excess in mandible.

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Fig. 17. Case with mandible Bolton Index "12" excess and no Bolton Index "6" excess.

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To sum up, it can be said that Bolton Index indicates where the excess of dental size is located and consequently where the stripping should be carried out, and where it should not be done. If there is no Bolton index discrepancy, the stripping should be done in both arches in order to not modify the relationship of size between the maxilla and mandible.

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Anatomic studies reveal that the minimal thickness of an enamel layer in the contact point is 1 mm in all teeth except upper lateral incisors and four lower incisors. As most of authors accept the rule which says that it is possible to remove half of the enamel layer thickness, it is acceptable to trim:

- 0.5 mm on mesial and distal surface of upper central incisors, upper and lower canines, upper and lower bicusps, upper and lower molars.
- 0.3 mm on mesial and distal surface of upper lateral incisors and four lower incisors.

	Central incisor		Lateral incisor		Canine		First bicuspid		Second bicuspid		First molar	
	mesial	distal	mesial	distal	mesial	distal	mesial	distal	mesial	distal	mesial	distal
Maxilla	0.5	0.5	0.3	0.3	0.5	0.5	0.5	0.5	0.5	0.5	0.5	0.5
Inter-dental reduction	1.0	0.8	0.8	0.8	1.0	1.0	1.0	1.0	1.0	1.0	1.0	1.0
Mandible	0.3	0.3	0.3	0.3	0.5	0.5	0.5	0.5	0.5	0.5	0.5	0.5
Inter-dental reduction	0.6	0.6	0.8	0.8	1.0	1.0	1.0	1.0	1.0	1.0	1.0	1.0

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Fig. 18. Chart with the recommended stripping values.

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Fig. 19. The Stripping Chart can also be used for the Stripping Treatment Plan.

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Facial and dental midlines

The facial midline is determined by glabella, tip of the nose, upper labial philtrum, Cupid's bow center and a chin. In an ideal esthetic situation, upper and lower dental midlines have to match between each other and with facial midlines.

The treatment goal is to achieve that the upper and lower dental midlines match with the facial midline.

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Fig. 20. Facial and dental midline.

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Teeth size

The teeth should be proportional among themselves and in relation to the patient's face. As Levin has demonstrated, the "Golden proportions" have always been a constant feature in the nature, art and design.

According to the "Golden proportion", the anterior teeth proportion from the frontal point of view should be:

- Upper Central incisor – 1.618
- Upper Lateral incisor – 1.0
- Upper Canine – 0.618


Fig. 21. Tamaño de dientes.

Fig. 21. Teeth size.

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But, if absolute measures are taken into account, Sterrett et al. determined the recommended measurements.

Tooth	Height (male)	Height (female)	Width (male)	Width (female)
Central incisor	10.2	9.4	8.6	8.1
Lateral incisor	8.7	7.8	6.6	6.1
Canine	10.1	8.9	7.6	7.1

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Fig. 23. Teeth shapes depending on age, sex and personality, according to Lombardi. The scheme shows the recommended teeth shapes depending on sex, age and personality, according to Lombardi.

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These measurements and proportions are useful when a decision has to be made on aesthetic retouches of the teeth, whether by stripping or by reconstruction.

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Design and planning of stripping

For the design and planning of the stripping it is necessary to take into account:

- Dento-alveolar discrepancy.
- Bolton Index "12" and "6".
- Position of the teeth in relation to the facial midline.
- Upper and lower dental midline discrepancies.
- Shapes of the teeth.
- Dental size proportions.
- Dental asymmetries.
- Intercuspidation.
- Distance between the contact point and the bone crest.
- If the teeth show a contact point or a contact surface.

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- Verde 90 μ
- Gris 60 μ
- Roja 40 μ
- Blanca 25 μ
- Amarilla 15 μ

Figs. 24 y 25. Echarri PST Set.

When the interproximal contact is a surface and not a point, it means that the teeth have already suffered an abrasion and this fact will limit the stripping possibilities.

The stripping should be carried out to meet all these goals if possible:

- Correct the crowding.
- Compensate the Bolton's discrepancies.
- Match the upper and lower dental midlines with the facial midline.
- Compensate the abnormal shapes of the teeth.
- Compensate the dental size proportions.
- Compensate the asymmetries of the teeth.
- Achieve the intercuspidation of upper and lower teeth
- Fix the inter-proximal contact point to 5 mm away from the bone crest to achieve the presence of the papilla
- Establish a correct surface of the interdental contact point.

Progressive Stripping Technique (PST)

The PST recommends using the files of the PST Kit to perform the stripping. It includes the following files:

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The upper central incisors, and all canines, bicuspsids and molars can be trimmed 0.5 mm and for this purpose we recommend to use:

File 60 μ during 60 seconds.
File 25 μ during 30 seconds.
File 15 μ during 30 seconds.


The upper lateral incisors, and all lower incisors can be trimmed 0.3 mm and for this purpose we recommend to use:

File 40 μ during 30 seconds.
File 25 μ during 30 seconds.
File 15 μ during 30 seconds.

The teeth with prosthetic crowns or interproximal fillings often need the reshaping of the interproximal contact point and the dental size by using:

File 90 μ as much as necessary.
File 25 μ during 30 seconds.
File 15 μ during 30 seconds.

These files also have other advantages. They can trim hard tissues without damaging the soft tissues and it is easier to limit the stripping to the contact point area using a file and not a bur.

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Fig. 26. Relationship of the stripping file with the contact point or contact surface. The figure shows how a flat file can trim only the most convex area of the mesial or distal dental surface, including the cases with an increased surface of contact point.

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
It is very important to position the files perpendicularly to the occlusal plane and parallel to the dental axis to make a correct trimming of the teeth.

The PST recommends the separation of the teeth before the stripping to have better access to the contact point area, as well as to keep the hand-piece in a fixed position allowing the file to make the stripping only by the movement of the hand-piece and always only under the irrigation. Avoiding the circular movement of the hand-piece, a better polishing of the dental surface will be achieved.

It is also very important to avoid the vertical oscillating movement of the files in order not to change the shape of the interdental space and the interdental sulcus. If the interdental space is increased, the possibility of gingival "Black triangles" will be also increased. Any modification of the interdental sulcus shape leads to the circulation of the food during the mastication. Avoiding the horizontal oscillating movement of the files, the shape of the labial and lingual spaces is preserved and the circulation of the food is not modified.

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Conclusions

The adequate design and planning of the stripping is very important to achieve all treatment goals from the esthetic and functional point of view. Besides a good planning of the stripping, it is also important to have the adequate instruments and to follow the clinical protocol.

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1. Aprile H, Figún ME Anatomía odontológica. Tercera edición. Editorial El Ateneo, Buenos Aires (Argentina) 1960.
2. Bennet JC, McLaughling RP Consideraciones sobre la forma de la corona de los incisivos en el tratamiento ortodóncico. Rev Esp Ortod 1997; (27):359-69.
3. Sheridan JJ The physiologic rationale for air-rotor stripping. J Clin Orthod 1997; (31): 609-12.
4. Tarnow D, Margner WS, Fletcher P The effect of the distance from the contact point to the crest of bone on the presence or absence of interproximal papilla. J Periodontol 1992; (63): 993-96.
5. Echarri P Clear Aligner. 1st edition. Ripano SA. Madrid, España 2013.
6. Jerrold L, Lowenstien LJ The midline: Diagnosis and treatment. Am J Orthod Dentofac Orthop 1990; (97): 453-62.
7. Johnston CD, Burden DJ, Stevenson MR The influence of dental to facial midline discrepancies on dental attractiveness ratings. Eur J Orthod 1999; (21): 517-22.
8. Levin EI Dental esthetics and golden proportion. J prosthet Dent 1978; (40): 244-52.
9. Ricketts RM The biologic significance of the divine proportion and Fibonacci series. Am J Orthod 1995; (1): 105-26.
10. Kokich VO, Kiyah HA, Shapiro PA Comparing the perceptions of dentists and lay people to altered dental esthetics. J Esthet Dent 1999; (11): 311-24.
11. Sarver DM, Ackerman MB Dynamic smile visualization and quantification and its impact on orthodontic diagnosis and treatment plan. In: Romano R, editor; Bichacho N. y Touati B associated editors. The art of smile. London (Great Britain): Quintessence Publishing; 2005. Pp 109-139.
12. Echarri P Ortodoncia en adultos. Enfoque actual. Revista Ripano 2012; 9(24): 18-24.
13. Echarri P Progressive stripping technique. (CE Course). Cleveeland, USA: Odontos World Media, 2005.
14. Echarri P Progressive Slenderizing Technique. En: Melsen B, editor. Adult Orthodontics. Oxford (Great Britain): Blackwell Publishing; 2012.

References